

LEARNER RESPONSE SYSTEM

The ultimate collaborative learning tool. Allows students to express themselves through words, numbers, symbols, Likert scales and more. Self-Paced Learning functionality encourages students to work at their own pace and level.

TECHNICAL SPECIFICATIONS

Range	100m (328ft) [Free Space]
LCD Interface	128 x 128 pixels, 36mm x 38mm (1.4in x 1.5in), 4 level grey scale, Reflective Polarizer
User Input	8 Intelligent 'hot keys', 4 'navigation keys' and 10 alphanumeric keys 3 keyboard layout options, Latin, Arabic or Cyrillic
Max Number of Devices per System	500
Batteries	3xAA Alkaline, replaceable via security screw
Battery Life	Typical usage 2 Years
Communication	Wirelessly with the ActivHub 2.4 GHz inserted in computer or ActivBoard 300 and 300 Pro
Languages	English, Dutch, French, German, Italian, Portuguese, Spanish, Finnish, Norwegian, Swedish, Danish, Turkish, Kazakh, Russian, Arabic, Brazilian Portuguese. Available languages dependant on variant.
Certifications	CE, FCC, RoHS
Overall Dimensions	73mm x 132mm x 31mm (2.9in x 5.2in x 1.2in)
Weight	205g (7.2oz)
Operating Temperature	0°C - 50°C (32°F- 122°F)
Storage Temperature	-10°C - 60°C (14°F - 140°F)
Pack Contents List	ActivExpression devices packed in a soft carrying case, ActivInspire Professional Edition DVD, 1 ActivHub, 1 security screwdriver, Quick Start Guide (batteries are included in each device)



GENERAL SPECIFICATIONS

Response Types	Multiple Choice, True/False, Y/N, Multiple Response, Put in Order, Gap Fill words, Likert scales, Confidence Checking & Alphanumeric response
Text Response	Words, phrases, numbers, symbols and punctuation (up to 120 characters)
Question Data Import	ExamView® XML importer
Question Data Export	Questions and results can be exported to Excel™
Working with other Applications	Interface (Express Poll/Wonderwheel) tool can “float” over most computer applications such as web browsers, Word™, Excel™
Question Preparation	Question Wizard allows for preparing questions with correct answers Questions can be timed Express Poll allows for on-the-fly questioning
Response Monitoring	Results can appear instantly or viewed at a later time Users can be allowed to modify their response Results can be captured either by user’s name or anonymously
Self-Paced Learning Questions	Questions are downloaded to and can be viewed on each device
Self-Paced Learning Question Preparation	Question Wizard includes settings for order, difficulty level and duration of questions
Self-Paced Learning Response Monitoring	Responses are immediately visible via a color-coded timeline

MINIMUM COMPUTER REQUIREMENTS

Operating Systems	Windows XP SP2 and above (XP/Vista/7) Mac OSX 10.4.11 – 10.6.1 (Tiger/Leopard/Snow Leopard) Linux Ubuntu 9.04, 9.10; Debian
PC	Pentium 4 – 1 GHz processor (800MHz for Vista), 512 MB of RAM, 1024x768 resolution, 1.5 GB of Free Disk space
Mac	Intel Processors (Universal Binary), 512 MB of RAM, 1024x768 resolution, 3.0 GB of Free Disk Space
Installation	Operation and use of ActivExpression devices requires ActivInspire software, included ActivHub and an available USB port on either the computer or ActivBoard. Importing ExamView® and XML files, Export questions and answers to Excel™ and Self-Paced Learning requires ActivInspire Professional Edition (included). User with administrator/privileged access rights required for software installation.

WARRANTY AND SUPPORT

ActivExpression	One-Year Warranty (Return-To-Base) 90-Day Telephone Support and One-Year Email Support KnowledgeBase: www.prometheankb.com
------------------------	--

ExamView® is a registered trademark of eInstruction, Inc.

Page 2

©2010 Promethean Ltd. All Rights Reserved. ActivExpression is a trademark of Promethean Ltd. Any third-party product and company names are mentioned for identification purposes only and may be trademarks of their respective owners. All dimensions +/- 5mm (0.2in). All weights +/- 0.9kg (2lb). www.prometheanworld.com Product Specifications are subject to change without notice. 12/09